

# SDS-Polyacrylamide Gel

## Electrophoresis

### (SDS-PAGE)

Determining the MW of Protein by SDS-PAGE



P<sub>38</sub>

# Objective



Aim :

Comprehend the principle and application of Sodium dodecyl sulfate Polyacrylamide Gel Electrophoresis (SDS-PAGE)

- Preparation of PAG needs a long time, so do it at first.

# Procedure



1. Reagent preparation:
2. Gel Cassette Assembly
  - (1) Clean and dry the glass plate, shelf etc.
  - (2) Assemble the gel apparatus
3. Preparation of the Gel and pour it into the sandwich layer of glass plate

# 3. Preparation of the Gel



(1) Prepare the Separating Gel (10%) 10 ml per gel

ddH<sub>2</sub>O 3.90 ml

30% Acrylamide stock solution 3.33 ml

1.5M Tris-Cl buffer (pH 8.8) 2.50 ml

10% SDS 100 μl

10% AP 100 μl

10% TEMED 100 μl

Mix and pour the gel using a pipet,  
allow to polymerize by overlaying gently with water.

# Classification of PAGE



- ❖ **Native PAGE**
- ❖ **Discontinuous PAGE**
- ❖ **SDS-PAGE**
- ❖ **IPG-IEF (Immobilized pH Gradient-IEF)**
- ❖ **2-Dimensional Electrophoresis**

# Principle of PAGE





1. The basis of electrophoresis (**vertical slab gels**):  
A charged protein will migrate toward the oppositely charged region in an electric field.
2. PAGE:  
Polyacrylamide: **Acrylamide** which forms a linear polymer, can be cross-linked with **N,N`-methylene bisacrylamide**, to form a gel matrix of **controlled pore size**. Polymerization is catalyzed by **free radicals**, generated by **ammonium persulfate (AP)** in the presence of **TEMED**.
3. Acrylamide stock solutions with varying ratios of **acrylamide and bis-acrylamide** are routinely used, to create different **pore size**.

# Formation of Polyacrylamide Gel




## Acrylamide






Polymerization

## Polyacrylamide Gel




## N,N'-methylene bis-acrylamide






**diluted solution**



**concentration —●— linker**



# The Structure of Polyacrylamide Gel



**gel —●— joint**



# Principle of PAGE




vertical slab gels



# Principle of SDS-PAGE

4. SDS-PAGE: the exact rate of movement of a particular protein depends on its size.
5. SDS coats proteins with an approximately uniform charge-to-mass ratio of (-) charge , and proteins are approximately uniformly shaped into spheres, thus, proteins separated by SDS-PAGE are denatured.
6. Disulfides in proteins are broken by addition of  $\beta$ -mercaptoethanol or 1, 4-dithiothreitol (DTT), so running proteins in SDS-PAGE are reduced condition.



## 7. A discontinuous gel system :

**The Stacking Gel:** approx. 10% of the volume of the total gel and a lower %(2.5-4.5%) acrylamide and a lower pH (6.8). charged molecules move freely and proteins in a sample should **accumulate** in stacks of closely spaced bands before encountering the separating gel.

**The Separating gel:** containing a higher % (7-15%) acrylamide and at a higher pH(8.8), proteins separate into discrete bands **based on size**.

The molecular weights of protein can be estimated by measuring the mobility of protein standards on the same gel.

# Buffer System



- ❖ **Buffer for Stacking Gel:** 0.5 mol/L Tris-Cl, pH6.8
- ❖ **Buffer for Separating Gel:** 1.5 mol/L Tris-Cl, pH8.8
- ❖ **PAGE running buffer (pH8.3) :**
  - 0.025 mol/L Tris-Cl
  - 0.2 mol/L Glycine

## 2X Loading buffer

100 mM Tris-Cl (pH6.8)

200 mM DTT

4% SDS

0.2% Bromophenol blue

20% glycerol

(DTT should be added just before the buffer is used, from 1M stocking solution)

# PAGE Procedure



1. Reagent preparation:

2. Gel Cassette Assembly

(1) Clean and dry the glass plate, shelf etc.

(2) Assemble the gel apparatus

3. Preparation of the Gel and pour it into the sandwich layer of glass plate

# 3. Preparation of the Gel



(1) Prepare the Separating Gel (10%) 10 ml per gel

ddH<sub>2</sub>O 3.90 ml

30% Acrylamide stock solution 3.33 ml

1.5M Tris-Cl buffer (pH 8.8) 2.50 ml

10% SDS 100 μl

10% AP 100 μl

10% TEMED 100 μl

Mix and pour the gel using a pipet,  
allow to polymerize by overlaying gently with water.

# 3. Preparation of the Gel



(2) Prepare the Stacking Gel	5.0%	5 ml per gel
ddH <sub>2</sub> O	3.4 ml	
30% Acrylamide stock solution	0.83 ml	
1.5M Tris-Cl buffer (pH 6.8)	0.63 ml	
10% SDS	50 μl	
10% AP	50 μl	
10% TEMED	50 μl	

decant the water overlayer, pour the stacking gel, Insert the comb and allow to polymerize.

# PAGE Procedure



(3) Once the gel has polymerized, the comb can be gently removed.

## 4. Sample preparation

- (1) Combine 20  $\mu$ l of sample (a unknown protein and a mix protein) with 20 $\mu$ l of protein loading buffer
- (2) Heat samples for 3-5 min in a boiling water bath to denature the proteins, then cold and **load the samples.**



# PAGE Procedure



## 5. Running the gel


(80V 0.5h; 150V 1.5V)

## 6. Coomassie Brilliant Blue Staining

(Staining and destaining)

## 7. Observation

(next week)



# PAGE Procedure




During running the PAGE,

Please see a video of SDS-PAGE.

# Result of SDS-PAGE



For example:



12% Separating Gel

# Discussion



- ❖ If given a sample containing two proteins with large different mass weight, please list the methods for separating two proteins, and briefly introduce the principles.