

武漢大学

Wuhan University

ELECTROPHORESIS



Dept. of Biochemistry Wuhan University

1001111











2

Supporting Mediums



Detection of Components after Electrophoresis (Dying)



Special Electrophoretic Techniques and Applications



Electrophoresis is the forced migration of charged particles such as macromolecules, in an electric field.
Cations move toward the cathode

and anions move toward the anode.



Biochemistry experiment 2





a voltage power supply, electrodes, buffer, and a supporting medium



To separate and purify the charged biomolecules:
 Protein, amino acids; Nucleic acid (DNA, RNA)



2016-4-7



For example: protein Electrophoresis





SEPARATION OF SERUM PROTEINS BY CELLULOSE ACETATE MEMBRANE (CAM) ELECTROPHORESIS



Master the principle and method of electrophoresis

How to separate and identify serum proteins using cellulose acetate membrane (CAM) as a supporting medium --CAM electrophoresis.



Separation of Serum Proteins by CAM Electrophoresis

Principle:

- In alkaline buffer solution, serum proteins are negatively charged. They will migrate toward anode in electric field.
- Because of the differences of net charge and molecular weight, the migration velocities of proteins are distinct.
- Proteins separated on film are fixed in position and stained with a proteinspecific dye, such as amido black. After staining and destaining, five zones indicating fractions including albumin, a_1 -globulin, a_2 -globulin, β globulin, and γ -globulin of serum proteins can be visualized.
- Their percentages can be estimated by scanning with densitometry.

Procedure

P32

- 1. Immerse a CAM strip into the electrophoretic buffer, take it out and absorb water with filter paper.
- 2. Add 3~5µl sample on the rough side of membrane at the site 2 cm from cathode, and put the membrane on the electrophoretic device with the rough side of membrane down, connect the electrode and the film with filter paper. Electrophoresis is carried out at a constant voltage of 80V for 50~60 min.



Principle of Electrophoresis





* f

Electrophoretic mobility, m
(an intrinsic property)
m = v / E = EQ/fE = Q / f

 Q (decided by Charge)
 Charge varies as a function of building block molecule composition and buffer pH .

based on Size (decided by molecular weight), Shape



1. Properties of the Molecules

Molecular Shape Molecular Size Net charge



A long, loose protein travels at a slower rate than a globular protein in a gel.



A smaller molecular has a faster U





- 2. Properties of the electrophoretic System
 Electric field strength: E
- Property of the support medium
 influence separation in 3 aspects:

 restrictions on mobility (adsorbing 吸附)
 electroendosmosis 电渗
 effect of molecular Sieve凝成的分子纬



3、 Electrophoretic buffer (缓冲余统) Solution pH $I = \frac{1}{2} \sum c_i z_i^2$ Ionic Strength Adherence of buffer

The most suitable Ionic Strength is about 20~200 mM.



Support Mediums

Paper or Starch

Cellulose Acetate Membrane

Agarose Gel

Polyacrylamide Gel (PAG)



1. Proteins

Special protein dye: Amido Black-10B; Coomassie brilliant Blue R-250 or G250; Silver staining;

fluorescence labelled (Cy2, Cy3 and Cy5) and Radiolabelled Proteins (C-14, H-3, P-32)

Immunochemical detection: Western Blotting



2. Enzymes

Appropriate enzymatic methods

3. Glycoproteins

Schiffs Reagent (Neutral Glycoprotein); Alcian blue (Acidic mucopolysaccharides)

4. Lipid

Sudan black, Oil red O

5. Nucleic acids

Ethidium Bromide or SYBR gold



Special Electrophoretic Techniques and Applications

* Isoelectric focusing (IEF)
* Capillary electrophoresis (CE)
* Two-dimensional polyacrylamide gel electrophoresis (2D-PAGE, or 2DE)

Pulsed-Field Gel Electrophoresis (PFGE)



SEPARATION OF SERUM PROTEINS BY CELLULOSE ACETATE MEMBRANE ELECTROPHORESIS



- Electrophoresis is carried out at a constant voltage of 80V for 50~60 min.
- After electrophoresis, take out membrane, immerse into the dye (Amido Black-10B) solution for 1 min. Then put it into destaining solution for rinsing.
- 4. Observe and judge the shade of color and width of the bands to estimate their quantity.

P32



Serum Protein Electrophoresis -- Normal Pattern



A: Pattern of separated serum proteins visualized in the CAM strip.
B: Densitometric scanning of the CAM strip converts bands to characteristic peak of 5 fractions.



Serum Protein

Proteins	Pl	MW
albumin 54~73%	4.8	69 kD
a ₁ – globulin 2.8~5.1%	5.06	200 kD
a ₂ – globulin 6.3~10.6%	5.06	300 kD
β – globulin 5.2~11%	5.12	90~150 kD
γ – globulin 12.5~20%	6.85~7.5	156~300 kD



- The concentration of total plasma proteins is about 60~80 g/L. The plasma proteins are a complicated mixture that includes not only simple proteins but also mixed or conjugated proteins (glycoproteins and lipoproteins).
- 2. The CAM electrophoresis can separate the plasma proteins into 5 fractions. Albumin is the major protein of human plasma (38~48 g/L, more than 50% of the total plasma proteins. The liver produces about 12 g of albumin per day.



- 3. Albumin has two main functions:
 - (1) maintaining colloid osmotic pressure of blood;
 - (2) transportation: acting as a carrier molecule for bilirubin, fatty acids, trace elements and many drugs.

4. As an indicator of liver and kidney disease, immune deficiencies, malignancies of the immune system, acute and chronic infection, genetic deficiencies, central nervous system disease and numerous other pathologies.





Figure 11–3. Abnormal protein patterns by agarose gel electrophoresis. Anode is to the left; fractions are as labeled in Figure 11–1; all samples are serum except for C. A. Inanition in an elderly patient with low total protein and markedly depressed albumin. B, Nephrotic syndrome with elevated α_2 -macroglobulin and beta lipoprotein. C, Urine from protein-losing nephropathy. D, Iron deficiency with elevated transferrin. E, Broadly elevated gamma and low albumin due to liver disease. F, Oligoclonal gamma fraction in patient with renal disease. G, Hypogammaglobulinemia.



Figure 11-4. Serum protein electrophoresis: clinicopathologic correlations. (Courtesy of Dr. A. F. Krieg.)